

VED REGISTRERING/SØK/FOLDER	DATA EXPLORER -> SØK OG EKSPORT			FOLDER MANAGER -> EXPORT TO EXCEL	
Feltnavn	Table	Field	Feltnavn i Excelfil ved eksport	Field	Feltnavn i Excelfil
Additive	Sample	Additive	C_ADDITIVE	C_ADDITIVE	Additive
Anatomical localization	Sample	Anatomical Loc	C_ANATOMICAL_LOC	C_ANATOMICAL_LOC	Anatomical Localization
Barcode	Sample	Barcode	LABEL_ID	LABEL_ID	Barcode
Biobank	Sample	Biobank	CLINICAL_TRIAL	CLINICAL_TRIAL	Biobank
Cold Ischemia Time	Sample	C Cold Ischemia Time	C_COLD_ISCHEMIA_TIME	C_COLD_ISCHEMIA_TIME	C Cold Ischemia Time
Collection Date	Sample	Collection Date	SAMPLED_DATE	SAMPLED_DATE	Collection Date
Collection Site	Sample	C Sampled Location	C_SAMPLED_LOCATION	C_SAMPLED_LOCATION	Sampled Location
Collection Time	Sample	C Collection Time	C_COLLECTION_TIME	C_COLLECTION_TIME	Collection Time
Donor Id	Sample	Donor Id	CLIN_TRL_SUBJECT_ID	CLIN_TRL_SUBJECT_ID	Donor Id
Fasting	Sample	C Fasting	C_FASTING	C_FASTING	Fasting
Flag Comment	Sample	Has Flags	HAS_FLAGS	HAS_FLAGS	Has Flags
Informed Consent (knyttet til prøven)	Sample	Informed Consent	T_IC (dokumentnr vises)	Informed Consent	Informed Consent
Original Sample	Sample	Original Sample	ORIGINAL_SAMPLE	ORIGINAL_SAMPLE	Original Sample
Parent Sample ID	Sample	Parent Sample	PARENT_SAMPLE	PARENT_SAMPLE	Parent Sample
Prep By	Sample	Prep By	PREP_BY	PREP_BY	Prep By
Processing Method	Sample	Processing Method	T_PROCESSING_METH	T_PROCESSING_METH	Processing Method
Project	Sample	Project	PROJECT	PROJECT	Project
Sample Amount	Sample	Sample Volume	SAMPLE_VOLUME	SAMPLE_VOLUME	Sample Volume
Sample Category	Sample	Sample Category	T_SAMPLE_CATEGORY	T_SAMPLE_CATEGORY	Sample Category
Sample Form	Sample	Sample Form	T_SAMPLE_FORM	T_SAMPLE_FORM	Sample Form
Sample ID	Sample	Sample ID	TEXT_ID	TEXT_ID	Sample ID
Sample Number	Sample	Sample Number	SAMPLE_NUMBER	SAMPLE_NUMBER	Sample Number
Sample State	Sample	Sample State	T_SAMPLE_STATE	T_SAMPLE_STATE	Sample State
Sample Type	Sample	Sample Type	SAMPLE_TYPE	SAMPLE_TYPE	Sample Type ID
Sampled By	Sample	Sampled By	SAMPLED_BY	SAMPLED_BY	Sampled By
Sampling Timepoint	Sample	C Sampling Timepoint	C_SAMPLING_TIMEPOINT	C_SAMPLING_TIMEPOINT	Sampling Timepoint
Snomed	Sample	C Snomed	C_SNOMED	C_SNOMED	C Snomed
SOP (benytt DE til søk*)	Sample	T Sop (ref no)	T_SOP	T_SOP (ref no)	SOP
SOP (benytt DE til søk*)	Sample	C Sop Reference (tittel)	C_SOP_REFERENCE	C_SOP_REFERENCE/SOP_USED (tittel)	C Sop Reference/Sop Used
SOP (benytt DE til søk*)	sample_view	Sop (ref no)	SOP		
SOP (benytt DE til søk*)	sample_view	Sop Used (tittel)	SOP_USED		
Thaw Count	Sample	T Thaw Count	T_THAW_COUNT	T_THAW_COUNT	Thaw Count
Time To Store	Sample	C Time to Store	C_TIME_TO_STORE	C_TIME_TO_STORE	Time to Store
Unit	Sample	Sample Units	SAMPLE_UNITS	SAMPLE_UNITS	Units ID
Visit	Sample	Visit	CLIN_TRL_VISIT	CLIN_TRL_VISIT	Visit
Warm Ischemia Time	Sample	C Warm Ischemia Time	C_WARM_ISCHEMIA_TIME	C_WARM_ISCHEMIA_TIME	C Warm Ischemia Time

*benytt begge SOP felt ved søk, for å sikre at man inkluderer all SOP-informasjon